


Rapid COVID-19 Diagnosis using Raman Spectroscopy and Machine Learning

Side Deck - Capstone 1



Author: Isaac
Ghebregziabher
Capstone Project One

July 4, 2021

Introduction

Current diagnostic
method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling

Summary

Acknowledgement

- ▶ > 171 million active cases
- ▶ 3.5+ million deaths
- ▶ Fast and reliable diagnostic is needed

SARS-CoV-2



>200 countries affected



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling


Summary

Acknowledgement

RT-PCR – Current COVID-19 detection method is time consuming and expensive

Side Deck - Capsone 1

- ▶ 3 days for sample preparation and RNA extraction
- ▶ Expensive PCR



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling


Summary

Acknowledgement

Principle of Raman effect

- ▶ Most light scatters unaffected (Rayleigh scattering)
- ▶ A few percent gets Raman scattered
- ▶ Raman Scattered light is signature of molecular composition

☐ Shine light on sample



- ☐ Most of the light is unaffected
- ☐ Small percentage of light undergoes frequency shift

Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis


Modeling

Summary

Acknowledgement

Rapid detection of Covid-19 using Raman spectroscopy and Machine Learning

Side Deck - Capstone 1



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling

Summary

Acknowledgement

Overview of Dataset obtained from Kaggle

Side Deck - Capsone 1

- ▶ 309 rows X 901 columns
- ▶ Each column → Raman Wavenumber
- ▶ Rows → Intensity for wavenumbers
- ▶ Each row corresponds to one observation
- ▶ Last column '*diagnostic*' is target variable

Feature names = Wavenumbers in units of cm-1

	400	402	405	407	410	412	415	417	420	422	...	diagnostic
0	0.0	-0.015237	-0.030607	-0.038309	-0.039078	-0.035809	-0.031176	-0.030395	-0.033311	-0.031603	...	Healthy
1	0.0	-0.012098	-0.028164	-0.035189	-0.036138	-0.031050	-0.026015	-0.027539	-0.028084	-0.027075	...	Healthy
2	0.0	-0.013000	-0.029058	-0.035021	-0.034994	-0.033025	-0.028413	-0.028470	-0.029737	-0.029198	...	Healthy
3	0.0	-0.015728	-0.034346	-0.045140	-0.047671	-0.044334	-0.040807	-0.040474	-0.041417	-0.040699	...	Healthy
4	0.0	-0.020355	-0.045839	-0.060556	-0.065805	-0.064988	-0.062097	-0.061955	-0.064759	-0.066886	...	Healthy

Raman scattered intensities for each Wavenumber

Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling

Summary

Acknowledgement

Raman detector might have dead pixels corresponding to all zero value intensities

Side Deck - Capsone 1

- ▶ No missing values in dataset
- ▶ 9 features wave-numbers with 0 intensity value
- ▶ Drop null features (treated as missing)




Figure: Single valued features (dropped).

Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling

Summary

Acknowledgement

- ▶ No class imbalance issues
- ▶ Dataset is balanced with $\approx 50:50$ class ratio.




Figure: $\approx 50 : 50$ COVID to Healthy class ratio.

Introduction

Current diagnostic
method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling


Summary

Acknowledgement

Visually indiscernible Raman spectrum

- ▶ Visually difficult to easily identify COVID from Healthy.
- ▶ Machine Learning model is needed for fast and reliable COVID diagnosis.

Average Spectrum: Healthy Vs CoVid



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling


Summary

Acknowledgement

High dimensional data visualization - PCA

Side Deck - Capsone 1

- ▶ Over 50 percent variance explained with two principal components.
- ▶ No visual class separation



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis


Modeling

Summary

Acknowledgement

Principal component analysis – feature reduction

- ▶ Over 90% data variance explained with 15 components.
- ▶ Feature reduction to 15 from 900!



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling


Exploratory Data Analysis

Modeling

Summary

Acknowledgement

Three models considered:



Introduction

Current diagnostic
method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling

Summary

Acknowledgement

Classification Report Training dataset: Logistic Regression

Misclassified samples: 6

Logistic regression: PCA 15 components

Classification report for train-dataset:

	precision	recall	f1-score	support
Healthy	0.98	0.96	0.97	112
SARS-CoV-2	0.97	0.98	0.97	119
accuracy			0.97	231
macro avg	0.97	0.97	0.97	231
weighted avg	0.97	0.97	0.97	231

97 percent accuracy!

Need to be tested with the test split.

Introduction

Current diagnostic
method

Raman spectroscopy
Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling

Summary

Acknowledgement

Classification Report Training split: Decision Tree

Misclassified samples: 0

Decision tree: PCA 15 components

Classification report for train-dataset:

	precision	recall	f1-score	support
Healthy	1.00	1.00	1.00	110
SARS-CoV-2	1.00	1.00	1.00	121
accuracy			1.00	231
macro avg	1.00	1.00	1.00	231
weighted avg	1.00	1.00	1.00	231

100 percent accuracy! Model might be over-fitting.

Need to be tested with the test split.

Introduction

Current diagnostic
method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling

Summary

Acknowledgement

Classification Report Training split: Random Forest

```
Misclassified samples: 0
Random Forest: PCA 15 components
Classification report for train-dataset:
              precision    recall  f1-score   support

   Healthy           1.00      1.00      1.00        110
  SARS-CoV-2         1.00      1.00      1.00        121

 accuracy                   1.00         231
 macro avg           1.00      1.00      1.00         231
 weighted avg        1.00      1.00      1.00         231
```

100 percent accuracy! Model might be over-fitting.

Need to be tested with the test split.

Introduction

Current diagnostic
method

Raman spectroscopy
Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling

Summary


Acknowledgement

Models Testing With Test Split

Side Deck - Capsone 1

- ▶ All models performed well
- ▶ RF and Logistic regression have highest accuracy

Confusion matrix: Test dataset



With PCA physical meaning of features is lost.

Need to model using all features without PCA to gain physical insight on features.

Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling


Summary

Acknowledgement

Random forest classifier is the best performing model

Though all models persisted with good accuracy:

- ▶ Random forest performs the best
- ▶ 97% classification accuracy
- ▶ Random forest chosen for deployment



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling


Summary

Acknowledgement

Important features with high predictive power

15 features out of 901 are the most important.

- ▶ Wavenumber in range [650, 870] has high predictive power
- ▶ feature 870 has the highest predictive power



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis


Modeling

Summary

Acknowledgement

Important features Raman band corresponds to RNA/DNA band

- ▶ Virus is an RNA/DNA protein
- ▶ Band [700, 900] is prominent for RNA/DNA
- ▶ Band corresponds to predicted important features



Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis


Modeling

Summary

Acknowledgement

Do we need more data to enhance model performance?

Model accuracy saturates well before the end of available data.



No need for more data.

- ▶ We developed supervised machine learning models COVID detection using Raman spectroscopy data.
- ▶ Logistic regression, decision tree, and random forest supervised machine learning algorithms were considered.
- ▶ We find COVID detection using Random forest results in highest detection accuracy of 97 percent.
- ▶ Future works needs to be done with covid suspect and covid survived data for more comprehensive and reliable conclusion.

Introduction

Current diagnostic method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data Analysis

Modeling

Summary

Acknowledgement

Springboard mentor: Yuxuan Xin

for time generous and insightful discussions

Introduction

Current diagnostic
method

Raman spectroscopy

Raman and ML

Data Wrangling

Exploratory Data
Analysis

Modeling

Summary

Acknowledgement